

AQUA OXY FASCIA THERAPY

■ HYDRATION ■ SKIN-OXYGENATION ■ ELASTICITY ■ FASCIA LIFT

Adopting the feature of air boxing rejuvenation strategy (skin oxygenation), and combine with EMS lifting current and RF deep thermal energy in one handle. Effectively solve the ageing problem from fascia layer and help to do collagen firming, muscle tightening, skin brightening and skin oxygenation at the same time.



2 in 1



- **Air Boxing** By encouraging the blood flow through quick and rhythmic air boxing, more oxygen and nutrition are available. It will achieve drainage and make the skin look radiant through the process of skin oxygenation.
- **EMS (electrical muscular stimulation)** can reduce fine lines and wrinkles by stimulating. Through EMS muscle lifting current to activate a rhythmic contraction of muscles, this helps muscle to do passive exercise and lifting.
- **RF Deep Thermal Energy** accelerates collagen fullness, skin tightening, wrinkle improvement, moisturize, elasticity restoration, tissue repair
- **Light Emitting Diode (LED) Therapy** Combined with the theory of LED therapy, red light could activate circulation and enforce the cell energy; blue light could produce a feeling of cooling and help the energy be balanced and concentrated.

Fascia AirGun




Focused AirGun



AQUA OXY FASCIA THERAPY


TREATMENT
CONCEPT A



Fascia AirGun

- Hydration Treatment
- Activation Treatment
- Pore Tightening
- Wrinkle Improvement

TREATMENT
CONCEPT B



Focused AirGun

- Fascia Lifting
- Elasticity Improvement
- Dark Circles / Puffiness

AQUA OXY FASCIA DEVICE



Single Function



2 in 1



3 in 1